•Please read carefully the entire instructions of this leaflet; it contains important information for the user.  
•Keep this leaflet; you may need to read it again.

ACC0057/1

**IALU** *Proct*<sup>®</sup>  
Hyaluronic Acid

### Form and presentation

IALU *Proct*<sup>®</sup> is in the form of a yellowish translucent greasy rectal ointment. It is available in a 30 g aluminum tube with a plastic rectal cannula.

### Composition

White petrolatum, paraffin Oil, propylene glycol, polidocanol, escin  $\beta$ -sitosterol phytosome, 18-  $\beta$  glycyrrhetic acid phytosome, butylated hydroxytoluene (BHT), hyaluronic acid sodium salt.

### Properties

IALU *Proct*<sup>®</sup> is a medical device (Class IIb) containing high molecular weight hyaluronic acid (in salt form), a fundamental component of extracellular tissue, which promotes healing by accelerating the process of repairing anal lesions.

Escin  $\beta$ -Sitosterol phytosome is a soothing, relieving and anti-itching botanical ingredient, which acts synergistically for effective symptom control.

18-  $\beta$  glycyrrhetic acid phytosome has decongestant properties and soothes the first symptoms of inflammation.

IALU *Proct*<sup>®</sup> ointment contains no steroids or topical analgesic ingredients, which makes it ideal for long-term use without any side effects.

### Indications

IALU *Proct*<sup>®</sup> ointment is indicated as an adjuvant treatment for the anal lesions' repair after proctological surgery, internal and external hemorrhoids, proctitis, cryptitis, anal rhagades, fissures and perianal fistulas.

### Mode and duration of use

For external use, apply locally a thin layer of ointment around the rectum with a clean finger.

For an internal application, a plastic rectal cannula is attached to the packaging:

1.Screw the provided cannula onto the tube, remove the protective cap and

squeeze the ointment until the cannula is full.

2.Insert the cannula gently and not too deeply into the rectum, and squeeze the tube gently as you slowly pull it out.

3.Stay still for 1 or 2 minutes

4.After each use, clean the outside of the cannula with a clean wet tissue or wet cotton and replace the protective cap.

It is recommended that you thoroughly wash your hands with soap and water before and after each use of this medical device.

Use IALU *Proct*<sup>®</sup> ointment twice a day, one in the morning and one in the evening, if possible after defecation.

### Contraindications

-Do not use IALU *Proct*<sup>®</sup> ointment in case of known hypersensitivity to any of its constituents. If local hypersensitivity reactions persist, consult your doctor.

### Interactions

Do not use in combination with other rectal medical devices or medicines without first consulting a doctor.

### Precautions for use / Warnings

- Prolonged use of any product for local use may cause sensitization. In this case, it is advisable to discontinue use of the treatment and consult a doctor to initiate appropriate treatment. If any other side effects occur, please inform your doctor or pharmacist.

- Keep out of the reach of children.

- Do not swallow.

- Avoid any contact with eyes.

- Close the product immediately after use.

- If you are pregnant or breastfeeding, ask your doctor or pharmacist for advice before using IALU *Proct*<sup>®</sup> ointment.

- The expiry date refers to a product in its closed packaging and stored under appropriate conditions. Do not use the product after the expiry date or if it is open or damaged.

- IALU *Proct*<sup>®</sup> is for strictly personal use.

- Do not dispose of the medical device in the sewage system or with household waste. Ask your pharmacist what to do with medical devices that are no longer used. These measures will help protect the environment.

### Storage conditions

Keep the medical device tightly closed away from light and heat, preferably at a temperature not exceeding 25 ° C.

**Text's revision date** : April 2025

Caution - Consult the precautions of use.

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Laboratoires Opalia Pharma S.A., Recordati group